English

Deutsch

Français

Español

Italiano

Nederlands

Português

User Manual **Weigh Modules** WXS/WXT

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1 Introduction

1.1 This User Manual

This User manual contains all information for the **operator** of the product.

- Read this User manual carefully before use.
- Keep this User manual for future reference.
- Pass this User manual to any further owner or user of the product.

1.2 Further documents

In addition to this printed User manual you can download the following documents from

Documentation WXS

▶ www.mt.com/ind-wxs-support

Documentation WXT

▶ www.mt.com/ind-wxt-support

- Technical data sheet
- Installation information (for trained personnel under the control of the operating company)
- Reference manual for the command set

1.3 Manufacturer information

The contact information of the manufacturer of the product is as follows:

- **Name:** Mettler-Toledo GmbH
- **Weblink:** <http://www.mt.com>
- **Physical address:** Im Langacher 44, 8606 Greifensee, Switzerland

2 Safety Information

2.1 Intended use

- Use the product only for weighing in accordance with this User manual.
- The weigh module is intended for indoor use only.
- Any other type of use and operation beyond the limits of technical specifications is considered as not intended.

2.2 Definition of signal words and warning symbols

Safety notes are indicated by signal words and warning symbols and contain warnings and information about safety issues. Ignoring safety notes can lead to personal injury, damage to the instrument, malfunctions and erroneous results.

Signal words

WARNING for a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.

CAUTION for a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data or minor or medium injuries if not avoided.

NOTICE (no symbol)
for important information about the product.

Note (no symbol)
for useful information about the product.

Warning symbols

General hazard

Electrical shock

2.3 Product specific safety notes

Your weigh module meets the state of the art technology and complies with all recognized safety rules, however, certain hazards can arise.

Do not open the weigh module: It does not contain any parts which can be maintained, repaired or replaced by the user. If you ever have problems with your weigh module, contact your authorized METTLER TOLEDO dealer or service representative.

Observe instructions

Always operate and use your weigh module only in accordance with the instructions contained in the product documentation. The instructions for setting up your weigh module must be strictly observed.

If the weigh module is not used according to the product manuals, protection of the weigh module may be impaired and METTLER TOLEDO assumes no liability.

Staff safety

Use only METTLER TOLEDO accessories and peripheral devices, these items are designed to work optimally with your weigh module.

Explosion hazard

It is not permitted to use the weigh module in explosive atmospheres of gases, steam, fog, dust and flammable dust (hazardous environments).

Safety notes

CAUTION

Risk of electric shock

The weigh modules may only be connected to DC power sources that meet the 12VDC +/-3% at all times.

The power supply must be approved by the respective national test center of the country in which the weigh module will be used.

3 WXS/WXT Weigh Modules

3.1 Specifications

Parameter		205	205DU	204
Maximum capacity		220 g	220 g	220 g
Readability		0.01 mg	0.1 mg	0.1 mg
Measurement properties				
Specification temperatures		10 ... 30 °C		
Specification humidity		20 ... 80% rH		
Limit values				
Maximum capacity, fine range		—	111 g	—
Readability, fine range		—	0.01 mg	—
Repeatability * (at nominal load)	sd	0.04 mg (200 g)	0.07 mg (200 g)	0.1 mg (200 g)
Repeatability, fine range * (at nominal load)	sd	—	0.03 mg (100 g)	—
Linearity deviation	sd	0.15 mg	0.2 mg	0.25 mg
Eccentricity deviation (test load)		0.3 mg (100 g)	0.3 mg (100 g)	0.4 mg (100 g)
Sensitivity offset ** (test load)		$2.5 \times 10^{-6} \bullet R_{nt}$	$3 \times 10^{-6} \bullet R_{nt}$	$4 \times 10^{-6} \bullet R_{nt}$
Sensitivity temperature drift		$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$	$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$	$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$
Sensitivity stability		$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$	$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$	$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$
Dynamics				
Settling time *** (typical)		3 s	2 s	2 s
Settling time, fine range *** (typical)		—	3 s	—

Parameter		26	26DU	3DU
Maximum capacity		22 g	22 g	3.2 g
Readability		0.001 mg	0.01 mg	0.01 mg
Measurement properties				
Specification temperatures			10 ... 30 °C	
Specification humidity			20 ... 80% rH	
Limit values				
Maximum capacity, fine range		—	11 g	1.2 g
Readability, fine range		—	0.001 mg	0.001 mg
Repeatability * (at nominal load)	sd	0.003 mg (20 g)	0.006 mg (20 g)	0.006 mg (3 g)
Repeatability, fine range * (at nominal load)	sd	—	0.0035 mg (10 g)	—
Linearity deviation	sd	0.02 mg	0.03 mg	0.02 mg
Eccentricity deviation (test load)		0.03 mg (10 g)	0.03 mg (10 g)	0.01 mg (2 g)
Sensitivity offset ** (test load)		$4 \times 10^{-6} \bullet R_{nt}$	$4 \times 10^{-6} \bullet R_{nt}$	$2 \times 10^{-5} \bullet R_{nt}$
Sensitivity temperature drift		$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$	$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$	$1.5 \times 10^{-6}/\text{°C} \bullet R_{nt}$
Sensitivity stability		$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$	$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$	$2.5 \times 10^{-6}/\text{a} \bullet R_{nt}$
Dynamics				
Settling time *** (typical)		3 s	3 s	3 s
Settling time, fine range *** (typical)		—	5 s	5 s

sd = Standard deviation
Rgr = Gross weight

Rnt = Net weight (sample weight)
a = Year (annum)

- * under good environmental conditions (free of vibration and air draft)
- ** applies only after adjustment at nominal capacity with an OIML E2 weight
- *** time between placing the weighed object on the weighing pan and indication of a stable weighing value under optimal environmental conditions

3.2 Unpacking the Weigh Module



Unpacking the electronic unit and terminal:

- 1 Remove the top foam cushion.
- 2 Remove the unit and the hardware.

Unpacking the load cell:

- 1 Remove the foam cushion and take out the hardware.
- 2 Carefully lift the white cushion containing the load cell out of the box and place it on a level surface.

- 3 Carefully pull the white cushioning material a few inches apart and remove the load cell.

Note

We recommend to keep the original packaging and to use it to transport the balance. The weighing pan must always be removed before transportation. Mount the protection cap to protect the interior of the weigh module.

3.3 Scope of delivery

Weigh module box

Parts	WXS/15	WXSS	WXTS	WXTP	WXTS3DU
WXS load cell (component version) with attached plastic cover to cover the weighing pan retainer.	✓	✓	—	—	—
WXT load cell (benchtop version) with attached plastic cover to cover the weighing pan retainer.	—	—	✓	✓	✓
Standard weighing pan	✓	✓	✓	✓	✓
Adapter weighing pan	✓	✓	✓	✓	—
Wind ring (only with WXS26 and WXT26 models)	✓	✓	✓	✓	✓
WXS/WXT User Manual (this document)	✓	✓	✓	✓	✓
Production certificate and CE Declaration of Conformity	✓	✓	✓	✓	✓

Weighing unit box

Parts	WXS/15	WXSS	WXTS	WXTP	WXTS3DU
WXS electronic unit (component version)	✓	✓	—	—	—
WXT electronic unit (benchtop version)	—	—	✓	✓	✓
Electronic unit-load cell connector cable, 0.5 m/1.6 ft long	—	—	✓	✓	✓
Electronic unit-load cell connector cable, 1.5 m/5 ft long	✓	✓	—	—	—
Mounting bracket for electronic unit, including clip and screws for attaching to DIN standardized rail	✓	✓	—	—	—
Terminal holder, including screws (for attaching the terminal to the electronic unit)	—	—	✓	✓	✓
AC adapter	✓	✓	✓	✓	✓
Power cable (country-specific)	✓	✓	✓	✓	✓

Terminal box

Parts	WXS/15	WXSS	WXTS	WXTP	WXTS3DU
SWT terminal (monochrome display), including protective covering	—	✓	✓	—	—
PWT terminal (color display), including protective covering	—	—	—	✓	✓
Terminal-electronic unit connector cable, 0.58 m/1.9 ft long	—	—	✓	✓	✓
Terminal-electronic unit connector cable, 2 m/6.5 ft long	—	✓	—	—	—
WXTS3 wind shield	—	—	—	—	✓

4 Operation

4.1 Electrical connections

WX weigh modules have an RS232C interface for the communication with the control systems.

Interface descriptions

Interface type	RS232C, bidirectional, full duplex, 600 to 38'400 bps
Interface update rate (max.)	23 weight values per second (with terminal) 92 weight values per second (without terminal)

Power supply for the weigh module

Input voltage	12 V DC nominal
Input current	0.84 ADC

- Use a stable power supply with no voltage fluctuations.
- If voltage fluctuations cannot be prevented, use a voltage regulator to deliver a constant voltage value to the load cell.
- The power supply must be approved by the respective national test center of the country in which the weigh module will be used.

Always operate the weigh module with the standard AC adapter delivered in the scope of delivery. The AC adapter is suitable for use with the following voltage range:

100 – 240 VAC, 50/60Hz

NOTICE

Check whether your local power supply falls within this range. If it does not, DO NOT connect the electronic unit or the AC adapter to your power supply and contact your local METTLER TOLEDO office.

4.2 Environmental conditions

WX weigh modules can be operated within the following environmental conditions:

Temperature Range	Operating / Compensated	+10 ... +30 °C
	Allowable ambient	+5 ... +40 °C
Relative air humidity		Max. 80 % at 31 °C, decreasing linearly to 50 % at 40 °C, non-condensing
Height above mean sea level		Max. 4'000 m (13'330 ft)
Warm-up time		At least 120 minutes after power up

4.3 Ingress protection

IP rating of the WX weigh modules:

When the plastic cover is attached and the connector cable is hooked up, the load cell achieves an IP45 rating.

If you wish to use a sprayer or water jet to clean the load cell, take appropriate measures to prevent water from building up between the base plate and the support. This could allow water to penetrate the load cell from below. Seal the gap between the support and the load cell base plate with a suitable gasket or sealing agent.

In operation, the load cell is rated IP30.

The electronic unit meets IP40.

The PWT and SWT terminals are rated IP54.

5 Maintenance

To ensure that your weighing module remains reliable, accurate, and functional for many years to come, the individual components must be cleaned and maintained periodically as appropriate for the intensity of use and the risk of contamination with debris.

5.1 Cleaning the weigh module

Clean the weighing pan and the load cell housing regularly with a damp cloth. The electronic unit and terminal can also be cleaned in this fashion as needed. For tougher dirt, a mild household cleaner may be used. Make sure that no liquid penetrates inside the components!

Please observe the following notes:

WARNING

Risk of electric shock

- 1 Disconnect the weigh module from the power supply prior to cleaning and maintenance.
- 2 Only use METTLER TOLEDO power cable, if these need to be replaced.
- 3 Ensure that no liquid comes into contact with the weigh module, terminal or AC adapter.
- 4 Do not open the weigh module, terminal or AC adapter. These contain no user-serviceable parts.


CAUTION

Damage to weigh module

Under no circumstances use cleaning agents containing solvents or abrasive agents, as this can damage or scratch certain surfaces (particularly the terminal screen).

You can also clean the load cell using a **low-pressure hand sprayer**. Before using a sprayer:

- 1 Remove the weighing pan and the wind ring and seal off the weighing pan retainer with the white plastic cover (1).
- 2 Make sure that the cover locks on properly!
- 3 Screw the connector cable to the electronic unit on.
- 4 Ensure that the jack (2) is sealed.
⇒ Under these conditions, the load cell meets IP45 requirements.
- 5 Dry the load cell it off with a soft cloth.

Note

The sealing function of the plastic cover may deteriorate over time with frequent use. Check the condition of the cover before you use it. If necessary, a new cover can be ordered as a replacement part **see Spare Parts** in the separate Installation Manual.

5.2 Cleaning the wind shield (WXTS3DU)

- 1 Open the left and the right wind shield door and remove the base plate protection from the wind shield.
- 2 Remove the top glass.
- 3 Use a soft cloth to clean the parts of the wind shield housing.

NOTICE

The elements of the wind shield can also be removed from the plate and be cleaned in a dish washer.

**CAUTION****Damage to the wind shield**

Do not use cleaning agents containing solvents or abrasive agents. This may damage the wind shield housing.

5.3 Maintenance

Your weigh module is a precision instrument, and periodic maintenance is one of the basic requirements to ensure it will perform well for you for many years to come.

Maintenance intervals will depend on the duration of use and the application and ambient conditions. Maintenance must be performed by a technician trained by METTLER TOLEDO.

Ask your METTLER TOLEDO office about service packages – regular maintenance by an authorized service technician will ensure your weigh module remains accurate and lengthen its service life.

5.4 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

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© Mettler-Toledo GmbH 08/2016